

PM-LIFE® OVO-PUNC

OOCYTE PUNCTURE NEEDLE

35



The PM-LIFE® OVO-PUNC oocyte puncture needle is used to perform ultrasound-guided puncture of the ovarian follicles for the purpose of an in vitro fertilisation.

ADVANTAGES


- Echogenic distal part for a good visualisation of needle positioning related to the follicle.
- Optimal puncture thanks to the length of the needle and the shape of its bevel.
- A mandrel obturates the needle during insertion to protect oocytes from any contamination

AVAILABLE REFERENCE(S)

Reference	Outer diameter (mm / Gauge)	Inner diameter (mm)	Useful length (cm)
PS1341530	1.5 / 17	1.2	30

ABOUT THIS MEDICAL DEVICE

1 st release date	2002
Single use	Yes
Class	Ila
Sterile	Yes
Sterilisation method	Sterilised via gamma radiation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 10 units in individual sterile wrapping
Storage	In its original box in a cool and dry place



PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES