

ENDOBIOPS


Endometrial sampling

Use

Allows performing hystologic sampling of the endometrium for cancer and hyperplasia testing

Endobiops Standard allows performing current sampling

Endobiops Plus allows improving the quality of the sampling



Endobiops Mini allows performing hystologic sampling for the patients showing a narrow uterine collar or cervix and on the post menopausal patients

Endobiops curve allows performing of laterally sampling for the patients showing a uterus with non-conventional form or to target the sampling.

Characteristics

Graduated transparent sheath, closed at the end and opening laterally holed at its distal part.

Piston allowing suction sampling.

Endobiops: Endometrial sampling

Reference	Designation	CH	Diam	Length	Packaging
PS1200923	Endobiops standard	9	3 mm	23 cm	25 per box
PS1200950	Endobiops LG50	9	3 mm	50 cm	10 per box
PS1210626	Endobiops Mini	6	2 mm	27 cm	25 per box
PS1220627	Endobiops Curved	6	2 mm	27 cm	25 per box

 **STERILE** **EO** Classe I sterile  0459